



DTS 90032955_001_03

K5 plus-Handstück 4911

Pos.	Mat.Nr.	Bezeichnung	Hinweis
	0.674.8172	Mutter	
	0.674.8222	Gehäuse	
	0.674.7392	Ring	
	0.411.3180	Pflegeset	
	0.224.7048	Tülle 5,5	
	0.070.3274	Schlauchleitung 3x0,5 LiF12Y11Y	
	0.674.8242	Druckstück	
	0.674.8472	KONTAKTGEHÄUSE 5,5	
	1.001.6751	Kappe K5 plus	
	1.001.6752	Motorgehäuse K5 plus	
	1.001.6740	Leitung K5 plus	
	1.001.6738	Griffhülse K5 plus	
	1.001.6737	Griffspitze K5 plus	
	0.223.0396	Kontaktbuchse	10 Stück / Pack
	0.257.4025	Kugel 4,0MM 3	
	0.223.0111	Kontaktstift	
	0.200.6250	O-Ring 15x1,5	
	0.200.6199	O-Ring 10,4x0,8	
	0.200.1045	Scheibe 10,5	
	0.200.0692	Filzscheibe 7x11x2	
	0.200.6360	O-Ring 18x2	
	0.200.6232	O-Ring 14x1	
	0.222.4164	V-RING JV16X1	
	0.201.7066	Druckfeder	
	0.220.0058	Kugellagereinheit 7x14,1x7,5	
	0.223.0082	Kabelstecker 4 pol	
	0.202.0021	Zylinderstift 2,35x40	
	0.220.0056	Kugellager 10x19x7	
1	0.674.4881	Spannzange 3/30	
1	0.674.4871	Spannzange 2,35/30	
2	0.674.8482	Abdeckung 3	

K5 plus-Handstück 4911

Pos.	Mat.Nr.	Bezeichnung	Hinweis
2	0.674.8212	Abdeckung	10 Stück / Pack
3	1.006.3347	A-Rotor 1.001.7669 K4+ / K5+	kein Vorab-Austausch möglich
3	1.001.7669	Rotor Ersatz	